

APPLICATION NOTE

Metals – Laser Flash Analysis



Measuring the Nearly Impossible – High-Precision LFA Measurements on Thin Copper: From Pulse Width to Signal Evaluation

Fabia Beckstein, Applications Laboratory, and Nico Dilsch, Software Department

Introduction

Only flash systems featuring high sensitivity, an appropriate pulse width, and advanced data evaluation can accurately measure thin, highly conductive materials. The biggest challenge when measuring such material is the extremely short measurement time. This requires both a high data acquisition rate and a very low pulse width.

Copper is a perfect example for this. With a thickness of 0.3 mm up to several millimeters, it is often used as a heat spreader, substrate layer or as a structured cooling plate, where both lateral heat distribution and reliable mechanical integration are required. Typical applications can be found in power electronics, battery technology and assemblies under high thermal stress, where compact design and efficient heat dissipation are crucial.

Method and Measurements Conditions

The LFA 707 *StratoFlash® Classic* is equipped with a laser that achieves high-energy density, which is particularly necessary at high temperatures. However, when measuring thin materials, low energy input is essential to prevent damage and overheating.

Thanks to its adjustable pulse width and voltage, the LFA 707 *StratoFlash® Classic* can adapt the energy input to the measurement requirements. The detector features a 2 MHz-data acquisition rate, ensuring a sufficient number of data points even at the shortest measurement times.

The measurement conditions are detailed in table 1.

Table 1 Measurement conditions

Material	Pure copper
Thickness	0.32 mm to 4 mm
Sample holder	Ø 12.7 mm
Temperature	Room temperature
Pulse width	100 to 600 µs
Model	Standard Model, based on Cape Lehmann with Pulse Correction

APPLICATIONNOTE Measuring the Nearly Impossible – High-Precision LFA Measurements on Thin Copper: From Pulse Width to Signal Evaluation

Measurement Results and Discussion


Figure 1 depicts the thermal diffusivity of copper of different thicknesses, ranging from 0.32 mm up to 4 mm. All results are within $\pm 2.5\%$ compared to the literature value of approx. $117 \text{ mm}^2/\text{s}$ at room temperature [1].

The pulse length was adjusted according to the thickness and measurement time, ranging from $100 \mu\text{s}$ to $600 \mu\text{s}$. The half time ($t_{1/2}$) varied over two orders of magnitudes from approx. $210 \mu\text{s}$ for the 0.32 mm sample to 24 ms for the thickest sample with 4 mm.


Figure 2 shows the signals for the samples of minimum and maximum thickness. The signal-to-noise ratio of both measurements is not ideal. This is due to the low energy input used to prevent overheating and the measurements being performed at room temperature. Nevertheless, the mathematical model fits the data perfectly, which is critical for achieving highly accurate

results. In laser flash analysis, the mathematical models used to determine the thermal diffusivity are based on the analytical solution of the heat conduction equation, assuming an instantaneous energy input (Dirac pulse). In reality, however, the laser pulse always has a finite duration. For samples with a relatively long measurement time, the pulse duration is typically much shorter than the characteristic measurement time, making deviations from the ideal assumption negligible (figure 2: 4 mm copper).

For highly conductive materials such as copper, especially when measuring thin samples, the thermal response occurs within a very short time. In such cases, the pulse duration is of the same order of magnitude as the characteristic diffusion time of the sample (figure 2: 0.32 mm copper). This leads to an overlap between the heating phase and the thermal response of the sample, which can distort the temperature curve and consequently, the calculated thermal diffusivity.




1 Thermal diffusivity of copper with different thicknesses at room temperature compared to literature values [1].



2 Detector signal (blue), mathematical fit (red) and pulse (orange) of 0.32-mm (left) and 4-mm (right) copper samples.

APPLICATIONNOTE Measuring the Nearly Impossible – High-Precision LFA Measurements on Thin Copper: From Pulse Width to Signal Evaluation



3 Thermal diffusivity of copper evaluated with and without pulse correction

Pulse Correction

To account for this effect, the NETZSCH LFA *Proteus*® analysis software automatically applies the exponential pulse correction [2]. Instead of assuming an instantaneous energy input, the real signal of the laser pulse is considered during the evaluation. This is achieved by incorporating the pulse signal through convolution, allowing the time-dependent heat input to be taken into account in the calculation of the temperature response. This way, the evaluated thermal diffusivity reflects the actual experimental conditions rather than an idealized instantaneous pulse.

By considering the actual pulse shape during evaluation, the pulse correction significantly improves the accuracy of the thermal diffusivity determination for thin and highly conductive samples. This becomes increasingly important as the sample thickness decreases and the thermal diffusivity increases.

For extremely short measurement times and thus also extremely short $t_{1/2}$, a robust and precise pulse correction is the most important analyzing feature. This is demonstrated in figure 3. As in figure 1, the blue dots represent the thermal diffusivity of copper with different thicknesses. In this case, the pulse correction was used for evaluation. The orange triangles represent the

same measurements, but the evaluation was performed **without** pulse correction. Decreasing the sample thickness – resulting in shorter measurement times – leads to increased errors caused by pulse overlapping.

Conclusion

The results demonstrate that even thin, highly conductive copper samples with extremely short thermal response times can be accurately measured using the LFA 707 *StratoFlash*® *Classic*. The combination of adjustable pulse control, high-speed data acquisition and advanced pulse correction ensures reliable thermal diffusivity results even under demanding measurement conditions. This makes the LFA 707 *StratoFlash*® *Classic* a powerful solution for the characterization of materials with very high thermal diffusivity.

Conclusion

[1] Touloukian, Y. S., et al. "Thermophysical properties of matter-The TPRC data series". Volume 10. Thermal Diffusivity. 1974

[2] Blumm, Opfermann, "Improvement of the mathematical modeling of flash measurements". High Temperatures-High Pressures, Volume 34(5):515-521, 2002