

Investigation of Cement Raw Material by STA-FT-IR

Dr. Carolin Fischer, Applications Laboratory Selb

Introduction

The thermal analysis of cement raw materials containing silicon dioxide, calcium carbonate, calcium sulfate dihydrate, and calcium hydroxide is a key approach to investigate the complex physical and chemical transformations that take place during heating and are decisive for clinker formation.

Simultaneous TGA-DSC measurements provide a combined view of mass changes and associated thermal effects, offering a comprehensive description of the material's overall thermal behavior across a wide temperature range. When complemented by FT-IR spectroscopy, the technique is further expanded by linking thermal events to the composition of the gases released

during heating, thereby significantly increasing the interpretative value of the analysis. In particular, the direct STA-FT-IR coupling based on the PERSEUS® concept provides clear advantages, as the FT-IR spectrometer is directly mounted on the STA furnace, resulting in a very short, heated gas path with minimal dead volume and excellent synchronization between thermal and spectroscopic signals, which is especially beneficial for the investigation of complex mineral systems. The small footprint of the coupled instrument setup fits into most lab surroundings.

Measurement Conditions

The measurement conditions are detailed in table 1.

Table 1 Measurement conditions

Instrument	NETZSCH STA <i>Jupiter</i> ® PERSEUS®
Temperature program	RT to 1450°C
Heating rate	20 K/min
Purge gas	Synthetic air, 70 ml/min
Crucible	Platinum, 85 µl, with lid and washer of Al ₂ O ₃ between the crucible and the sensor
Sample mass	24 mg

Results and Discussion

In the TGA-DSC diagram shown in figure 1, a sequence of several thermal processes can be identified that are typical for cement and cement related raw material and extend over the entire temperature range up to approximately 1400°C.


In the temperature range between 100 and 200°C, a mass loss of about 7.5% is observed in the TGA signal, accompanied by a DTG minimum at 149°C and two overlapping endothermic DSC effects with peaks at 153°C and 168°C. This region is characteristic for the release of physically bound water as well as the dehydration of calcium sulfate dihydrate to hemihydrate and/or anhydrite.

Between 400°C and 600°C, a further mass loss of approximately 3.5% occurs, associated with a DTG signal at around 453°C and an endothermic DSC peak with a peak temperature of 457°C. This behavior is typical for the dehydroxylation of calcium hydroxide, during which structurally bound water is released.

The effect observed in the DSC signal at approximately 575°C is characteristic for the reversible α - β phase transformation of quartz (SiO_2).

Between 700°C and 850°C, an additional mass loss of 5.9 % is detected, correlating with a clear DTG minimum at 779°C and an endothermic DSC signal with a peak temperature of 784°C. This step is characteristic for the thermal decomposition of calcium carbonate, i.e., decarbonation accompanied by the release of CO_2 .

The DSC Effect at 1216°C is a hint on a phase transition, marking the formation of silicate phases.




1 Temperature-dependent mass change (TGA, green), rate of mass change (DTG, black) and heat-flow curve (DSC, blue) of cement raw material.

APPLICATIONNOTE Investigation of Cement Raw Material by STA-FT-IR


Above approximately 1250°C, a mass loss of about 17% is observed, accompanied by several intense DSC signals with maxima at around 1318°C and 1386°C, as well as pronounced DTG peaks at 1321°C and 1386°C. Among other processes, the decomposition of CaSO_4 to CaO and the associated release of sulfur oxides occur in this temperature range. In addition, these effects also mark the transition from pure decomposition reactions to high temperature induced phase transformations and the onset of melting processes, which are typical for cement and clinker related system

The complete IR data is shown in figure 2 in temperature- and wavenumber-dependent 3D plot. The TGA curve is plotted in red at the back and shows the correlation of the mass loss to the increase in IR intensity. For detailed evaluation of the IR data, single IR spectra were taken at different temperatures and compared to the EPA-NIST library.



2 3D plot of all detected IR spectra of cement raw material, TGA curve plotted in red at the back of the cube

APPLICATIONNOTE Investigation of Cement Raw Material by STA-FT-IR



3 Temperature-dependent mass change (TGA, green) and traces of H₂O (turquoise), CO₂ (black) and SO₂ (red) of cement raw material.

This revealed the release of water during the first two mass-loss steps, which correlates well with the dehydration of calcium sulphate and dehydroxylation of calcium hydroxide. The release of carbon dioxide was found between 550°C and 800°C due to the decomposition of carbonates. The last mass-loss step released SO₂ from the sulphate decomposition. The traces of gas release can be easily correlated to the TGA curve, see figure 3.

Summary

The STA-FT-IR analysis of cement and cement related raw material enables a comprehensive characterization of the physical and chemical processes occurring during heating. By combining TGA and DSC, mass changes

and the associated thermal effects are recorded simultaneously, while FT-IR coupling allows the unambiguous identification of the gases released during these processes. This makes it possible to clearly assign individual reaction steps such as dehydration, dehydroxylation, decarbonation, and sulfate decomposition. A key advantage of the method is the direct correlation between mass loss, thermal effects, and gas composition, which significantly reduces ambiguity in the interpretation of overlapping reactions.

STA-FT-IR therefore represents a powerful tool for the analysis and optimization of cement raw materials and clinker formation processes.