



KINETICS
SOFTWARE



NETZSCH



Thermokinetic Analysis

NETZSCH reaction kinetics software permits the simultaneous evaluation of up to five dynamic test data sets. For the first time a reliable estimation of kinetic parameters and of complex reaction models


(multiple-step reaction) is possible. Incorporating advanced analytical techniques and statistical criteria, it represents the state-of-art in kinetic analysis.

NETZSCH reaction kinetics software permits the analysis of complex reactions such as competitive, parallel, independent and single step type reactions.


single-step reactions



double-step reactions



triple-step reactions




NETZSCH reaction kinetics software is:

- versatile** evaluates thermogravimetric, differential calorimetric, dilatometric, evolved gas data
- discerning** determines kinetic parameters of complex reactions
 - single step
 - double step
 - triple step
 with following, competitive, parallel and independent reactions
- reaction types**
 - n-th order
 - phase boundary
 - nucleation
 - diffusion
 - autocatalytic

- precise** utilizes non-linear regression, time-based 4-th order integration, correlation analysis, f-test for selection of the most probable kinetic model
- fast** offers pull-down menu driven format, direct access to all NETZSCH data sources by means of ASCII-files
- encompassing** provides complete graphical presentation
 - graphic of multiple-step reactions
 - graphic of FRIEDMAN- or OFW-Analysis
 - simulated engineering graphics

Model-free Estimation of Activation Energies:


FRIEDMAN-Analysis




The FRIEDMAN-analysis indicates that the pyrolysis of poly (1.4-phenylene-di(methoxy-phenyl) vinylene) in Ar flow proceeds as three-step process with activation energies of approx. 200, 263, and 258 kJ/mol for the 1st, 2nd and 3rd step, respectively. Heating rates: 1.0, 2.5, 5.0 and 10 K/min

For $\alpha = \text{const}$ the plot of $\lg \frac{d\alpha}{dt} = \lg A - 0.434 \cdot E/RT + f(\alpha)$ with the slope $m = -0.434 \cdot E/RT$. The values of $\lg(A \cdot s)$ are estimated assuming $f(\alpha) = 1 - \alpha$.

Alpha	Activation Energies		E kJ/mol	lg (A·s)
0.95	272.79	+ -	22.62	13.77
0.90	257.67	+ -	7.39	13.28
0.80	266.40	+ -	7.16	14.31
0.60	295.91	+ -	16.66	17.56
0.40	265.35	+ -	4.97	16.01
0.30	261.46	+ -	4.86	16.14
0.20	318.03	+ -	23.21	22.06
0.10	246.19	+ -	5.46	16.80
0.05	201.98	+ -	4.65	13.13
0.02	187.04	+ -	7.08	11.80



OZAWA-FLYNN-WALL-Analysis



OZAWA-FLYNN-WALL-analysis of the pyrolysis of poly (1.4-phenylene-di(phenoxy-phenyl)vinylene) in Ar Heating rates: 1.0, 2.5, 5.0 and 10 K/min


In experiments with constant but different heating rates, $\beta = dT/dt$, the plot of $\lg \beta$ vs. $1/T$ of points with the same degree of conversion $\alpha = \alpha$ yields lines with the slope $m = -b \cdot E/RT$. The value of $b = 0.457$, given by DOYLE, is correct for $E/RT \approx 37$. Hence, for a more precise estimation of activation energies the NETZSCH kinetics software considers the dependence of b from E/RT by a new straightforward strategy.

Engineering Graphics

The engineering graphics allow predictions on the isothermal behaviour of the sample investigated by non-isothermal methods.


In this way advantageous conditions for the curing of a phenolic resin can be found, for example.

Isothermal plot of degree of conversion vs. time for the curing reaction of a phenolic resin. For this example a 1st order reaction with autocatalysis is obtained in the kinetic analysis.



parameters: $\lg 10 A1\text{-s}$: 6.70
 E/kJmol^{-1} : 77.86
 $\lg k\text{-cat}$: 0.27


Isoconversion plot for different degrees of conversion vs. time and temperature shows the conditions to achieve a certain progress of the reaction beside 100% conversion. Here the conditions for the vulcanization of an EPDM rubber are given according to an n-th order reaction model.



parameters: $\lg 10 (A1\text{-s})$: 23.00614
 $E1 / \text{kJmol}^{-1}$: 220.96712
 React. ord.1 : 1.19089

Isothermal plot of reactands vs. time, using kinetic parameters and a model obtained by non-linear regression of the pyrolytic decomposition of poly(1,4-phenylene-di(phenoxy-phenyl) vinylene).

Temperature: 480°C




Analysis of multiple-step reactions by non-linear regression (NLR)

The precise curve fitting of data obtained under different conditions is the reliable base

for each prediction. The advantages of the multiple-scan analysis in combination with

NLR are both high flexibility of model selection and relevance of the estimated results.


The pyrolytic decomposition of a poly(1,4-phenylene-di(phenoxy-phenyl)vinylene) sample in Ar-flow is identified as a follow-up reaction when evaluating TG mass loss curves. Heating rates: 1.0, 2.5, 5.0 and 10 K/min



Parameters and deviations

Parameter	Optimum value	t-Stand. Dev.
$\lg 10 (A1\text{-s})$	11.38398	0.05819
$E1 \text{ kJ/mol}$	228.62285	0.95264
$\lg k\text{-cat } 1$	1.08848	0.01865
$\lg 10 (A2\text{-s})$	15.28171	0.50598
$E2 \text{ kJ/mol}$	282.99333	4.17640
React. ord. 2	2.20523	0.04647
Foll. React. 1	0.51529	0.04774
Mass Diff. 1/%	-100.00000	const.
Mass Diff. 2/%	-100.00000	const.
Mass Diff. 3/%	-100.00000	const.
Mass Diff. 4/%	-100.00000	const.

The vulcanization of an EPDM rubber was measured by DSC at five different heating rates (1 to 20 K/min). The exothermal DSC peaks are recorded from 130°C to 220°C. The curve fitting in the multiple scan analysis by the non-linear regression method (NLR) is excellent. The kinetic parameters show an overlapping endothermal effect (negative following reaction) upon the main curing reaction.



Parameters and deviations

Parameter	Optimum value	t-Stand. Dev.
$\lg 10 (A1\text{-s})$	18.57249	0.03046
$E1 \text{ kJ/mol}$	178.18989	0.41617
React. ord. 1	0.98844	0.05809
$\lg 10 (A2\text{-s})$	0.75647	0.25688
$E2 \text{ kJ/mol}$	19.76593	2.36022
React. ord. 2	0.77228	0.02542
Foll React 1	-0.32654	0.03605
Area1/(J/g)	15.05956	const.
Area2/(J/g)	16.59976	const.
Area3/(J/g)	15.45947	const.
Area4/(J/g)	19.46043	const.
Area5/(J/g)	15.11410	const.

Fields of Application

Thermogravimetry (TG)

decomposition
oxidation
adsorption
desorption

Differential Scanning Calorimetry (DSC)

decomposition
oxidation
polymerization
curing/cross linking
vulcanization
crystallization

Thermomechanical Analysis/ Dilatometry (TMA/DIL)

sintering

subject to technical change

Service of the Software

- prediction of technical processes (process optimization by fixing of time-temperature correlations and temperature cycling)
- test modelling
- stability and durability prediction
- information about intermediate-products of multiple-step reactions

Hardware-Requirements

Minimum system: IBM® AT-compatible with co-processor; 640 KB RAM; 2 MB harddisk space; MS DOS version 3.1 or newer; CGA, EGA, Hercules or VGA graphics adapter.

Recommended system: 386 SX/386 with co-processor or 486; Epson graphics compatible 8/9 or 24 pin oder ink jet printer.

Further Software Solutions

Thermal analysis software for the acquisition and analysis of TG, DSC, TMA/DIL data. Thermal analysis mass spectrometer software for the programming, data controlling and evaluation of simultaneous TA-mass spectrometer data.

**NETZSCH –
The Exact Solution**

NETZSCH
Thermal Analysis

Headquarters:
NETZSCH-Gerätebau GmbH
D-95088 Selb/Bavaria
P. O. Box 1460
Wittelsbacherstraße 42
Phone: (0 92 87) 8 81-0
Fax: (0 92 87) 8 81 44

NETZSCH do Brasil
Industria e Comercio Ltda.
Rua Michigan, 166 Brooklin
BR-04.566 000 Sao Paulo
Tel.: (011) 536 9166
Fax: (011) 543 2716

NETZSCH Italiana S.r.l.
Sezione Strumentazione
Via Pasteur 10/12
I-37135 Verona
Tel.: (045) 820 1462
Fax: (045) 820 1463

NETZSCH Instruments, Inc.
4 Landmark Square
Stamford, CT 06901, USA
Tel.: (203) 358 8686
Fax: (203) 358 9142

NETZSCH-Gerätebau GmbH
Representative Office Singapore
106, Tagore Lane Block 3
Singapore 2678
Tel.: (065) 457 2103
Fax: (065) 451 6894

NETZSCH-Gerätebau GmbH
Bureau de Liaison de Lyon
L'Orée d'Ecully, Allée A
Chemin de la Forestière, F-69130 Ecully
Tel.: 78 33 13 53
Fax: 78 33 48 90

NETZSCH-Gerätebau GmbH
NETZSCH Mastermix Ltd.
Vigo Place, Aldridge, Walsall
West Midlands WS9 8UG, England
Tel.: (1922) 590 06
Fax: (1922) 533 20